

Rail clamping system · Wall brackets · Stainless steel rail

Rail clamping system



Universal bracket made of anodized aluminium for stand Ø 30 mm.



ABS clamp for rail 30x10 mm and 50x10 mm complete with nipple threaded ISO G. 1/4" F. and hose connector Ø 6 mm (for medical gases) or Ø 7.5 mm (for Vacuum).



ABS clamp bracket for rail 30x10 mm and 50x10 mm, with slide 25x5 mm, 30x5 mm or 45x5 mm.



Chrome-plated universal bracket for rail 30x10 mm with locking knob.



Chrome-plated universal bracket for rail 30x10 mm and 50x10 mm.



Anodized aluminum wedge for OHMEDA rail clamp bracket, with nipple threaded ISO G. 1/4" F. and hose connector Ø 6 mm (for medical gases) or Ø 7.5 mm (for vacuum).



Anodized aluminum OHMEDA rail clamp bracket.



Anodized aluminum wedge for OHMEDA rail clamp bracket, with slide 25x5 mm, 30x5 mm or 45x5 mm.

< follows

Rail clamping system · Wall brackets · Stainless steel rail

Wall brackets



ABS wall bracket: 25x5 mm, 30x5 mm, 45x5 mm slide.

Stainless steel rail



Directory, components detail

- 1 Chemical fixings
- 2 Plasterboard fixings
- 3 Rail spacer
- 4 Rail
- 5 Bolt washer
- 6 Bolt

Stainless steel rail 30x10 mm,
brushed, with spacers and fittings,
various lengths.

