

SUPPORTING DEVICES

FOR COLLECTION JARS

The range of jars supporting devices includes trolleys, metallic floor stands, ABS support rings and support baskets for all containers sizes and models.

Baskets and rings are manufactured with a strong and stiff hook matching with main slide sizes.



ABS SUPPORT RING FOR **MAK/1000** JAR,
SLIDES 25x5, 30x5, 41x4,
45x5 mm HOOK



ABS SUPPORT RING FOR **MONOKIT®**
AND **MAK/2000** JARS, SLIDES 25x5, 30x5,
41x4, 45x5 mm HOOK



ABS SUPPORT RING FOR **FLOVAC®** JARS,
SLIDES 25x5, 30x5, 41x4,
45x5 mm HOOK



METALLIC SUPPORT BASKET
FOR **MAK/4000** JAR, SLIDE 25x5 mm



CLIP FIXING 25x5 mm RINGS FOR TROLLEYS,
BRACKETS AND WALL PLATES



FLOOR STAND
SINGLE BASKET FOR **MONOKIT®**
AND **MAK/2000** JARS



5 WHEELS TROLLEY UP TO
4 DIFFERENT JARS
(MAK, MONOKIT®, FLOVAC®)
WITH VACUUM GAUGE AND 4 ON-OFF TAPS



5 WHEELS TROLLEY UP TO 4 DIFFERENT JARS
(MAK, MONOKIT®, FLOVAC®) WITH HANDLE



5 WHEELS TROLLEY UP TO 4 DIFFERENT JARS
(MAK, MONOKIT®, FLOVAC®) WITH VACUUM
REGULATOR AND 4 ON-OFF TAPS



5 WHEELS TROLLEY UP TO 4 DIFFERENT JARS
(MAK AND MONOKIT®) AND UP TO 12 JARS
(FLOVAC®) WITH VACUUM REGULATOR AND
SAFETY JAR EASYSAFE® AND 4 ON-OFF TAPS