

## Probes with thread connection



**AFNOR NF-S 90-116**  
probe, thread  
ISO G. 1/4" F.



**AFNOR NF-S 90-116**  
type **EASYFIX**® probe,  
thread ISO G. 1/4" F.



**UNI 9507** probe,  
thread ISO G. 1/4" F.



**BS 5682** probe,  
thread ISO G. 1/4" F.



**DIN 13260-2** probe,  
thread ISO G. 1/4" F.



**SS 875 24 30** probe,  
thread ISO G. 1/8" M.  
or ISO G. 1/4" F.



**JIS T 7101** probe,  
thread ISO G. 1/8" M.



**SANS 1409** probe,  
thread ISO G. 1/4" F.



**OHMEDA** probe,  
thread ISO G. 1/4" F.


# Probes with hose connection

## Flowmeter outlet hose connection · Stainless steel rail

### Probes with hose connection



### Flowmeter outlet hose connection



### Stainless steel rail




# Rail clamping systems for oxygen flowmeters

## Cylinder inlet connections for oxygen pressure regulators


### Rail clamping systems



ABS clamp bracket for rail 30x10 mm and 50x10 mm, with nipple threaded ISO G. 1/4" F. and hose connector Ø 6 mm.



Anodized aluminum wedge for OHMEDA rail clamp bracket, with nipple threaded ISO G. 1/4" F. and hose connector Ø 6 mm.



Anodized aluminum OHMEDA rail clamp bracket.

### Cylinder inlet connections for oxygen pressure regulators



EN 850 (or CGA 870) "pin-index" yoke cylinder connection.



AFNOR NF-E 29-650 cylinder connection.



DIN 477-1 cylinder connection.



UNI 11144 cylinder connection.



BS 341-3 "bull nose" cylinder connection.



CGA cylinder connection.

