

## Rail clamping system · Wall brackets · Stainless steel rail

### Rail clamping system



Universal bracket made of anodised aluminium for stand Ø 30 mm.



ABS clamp for rail 30x10 mm and 50x10 mm complete with nipple threaded ISO G. 1/4" F. and hose connector Ø 6 mm (for medical gases) or Ø 7.5 mm (for Vacuum).



ABS clamp bracket for rail 30x10 mm and 50x10 mm, with slide 25x5 mm, 30x5 mm or 45x5 mm.



Chrome-plated universal bracket for rail 30x10 mm with locking knob.



Chrome-plated universal bracket for rail 30x10 mm and 50x10 mm.



Anodized aluminum wedge for OHMEDA rail clamp bracket, with nipple threaded ISO G. 1/4" F. and hose connector Ø 6 mm (for medical gases) or Ø 7.5 mm (for vacuum).



Anodized aluminum OHMEDA rail clamp bracket.



Anodized aluminum wedge for OHMEDA rail clamp bracket, with slide 25x5 mm, 30x5 mm or 45x5 mm.

< follows

## Rail clamping system · Wall brackets · Stainless steel rail

### Wall brackets



### Stainless steel rail

