

# Accessories



## Jointed extensions

**flow-meter™** jointed extension is a system particularly suitable to support the corrugated tubes of respiratory circuits.

The jointed extension is available with or without rail clamp bracket, with one or two joints allowing an easy orientation and with a device to fix two tubes Ø 16 mm or Ø 22 mm. Thanks to the swivel joints with cogging, it is an efficient support for even considerable weights.

### Sizes of sections | Type with:

<b>ONE JOINT</b>	First section: 150 mm - Second section: 450 mm
<b>TWO JOINTS</b>	First section: 150 mm - Second section: 450 mm Third section: 300 mm



## Phleboclysis rods

**flow-meter™** phleboclysis rods are support units particularly suitable to hold the infusion bottles during the hospital therapies. The phleboclysis rods allow an easy orientation and they can have the two following constructions:

- fixed construction, with straight or angled terminal piece and two hooks;
- extendible construction, with angled terminal piece and four hooks.

Thanks to the particularly strong structure made of chrome-plated steel, they can support even considerable weights.

Different constructions are produced under request.